Plasmid spin miniprep kit manual (up to 30 ug)

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**Lumiprobe Corporation**  
(US and Worldwide)  
9:00AM - 9:00PM EST

201 International Circle, Suite 135  
Hunt Valley, Maryland 21030  
USA

Phone: +1 888 973 63 53

**Lumiprobe GmbH**  
(Germany and Europe)  
9:00AM - 5:00PM CE(S)T

Feodor-Lynen-Strasse 23  
30625 Hannover  
Germany

Phone: +49 511 16596811

**Biotech Industry Ltd**  
(Russia and CIS)  
9:00 - 18:00 GMT+3 (Moscow)

Leninskie Gory, 1, bld. 75  
119992 Moscow  
Russian Federation

Phone: +7 495 973 6353

ISO 9001:2015

The certificate No.: 12 100  
52989 TMS  
Issued by TÜV SÜD  
Management Service GmbH

# Plasmid spin miniprep kit manual (up to 30 ug)

## Plasmid DNA purification protocol

### Before you start:

Add RNase A to the resuspension solution and mix (for more complete transfer of RNase A, wash the RNase A tube with a small amount of resuspension solution). To make wash solution concentrate ready add 40 ml of ethanol to 10 ml of wash solution concentrate and mix. Put marks on the lids of vials with a resuspension and washing solution.

Make sure that there is no precipitate in the lysis and neutralising solutions. If precipitate is visible, heat to 50°C in a water bath until the precipitate is completely dissolved.

Do not use more than 30 optical units of cultures for plasmid isolation (eg 15 ml a culture with OD600 equal to 2). If necessary to use more cells for plasmid purification you should proportionately increase the volumes of resuspension, lysis and neutralising solutions, and apply the supernatant several times on the column.

1. Centrifuge 2.5-7 ml of an overnight culture of bacteria (10-15 ml for low copy plasmids) at 5,000 rpm (4500g) for 5 minutes (or at 13000 rpm (<10000g) for 1 minute).
2. Carefully remove the supernatant. Resuspended the pellet in 250  $\mu$ l of resuspension solution and transfer it to a 1.5 ml tube.
3. Add 250  $\mu$ l of the lysis solution to the cell suspension, mix carefully by inverting the tube about 6 times (In case of using 10-15 ml of the culture, the five inverting may not be enough. In this case, the contents of the test tubes must be turned several times rapidly. Do not use vortex for mixing. The solution after this step should not be cloudy, it must be translucent and viscous.)
4. Add 350  $\mu$ l of neutralizing solution, gently mix by inverting the tube 10 times (if using 10-15 ml of the culture, shake the tubes for a few seconds; white flakes are formed in solution after mixing).
5. Centrifuge 10 000 - 13 000 rpm for 5 minutes.
6. Place the columns into the collection tubes. Transfer the supernatant into the column.
7. Centrifuge 10 000 - 13 000 rpm for 30 to 45 seconds. Discard the flow-through and place the column back into the same collection tube.

*Optional. To completely remove the nuclease, when isolating from end A + strains and removing endotoxins, add 500  $\mu$ l washing solution with guanidium chloride (use of this solution may result in a 20% reduction in yield). Centrifuge 10 000 - 13 000 rpm for 60 seconds. Discard flow-through.*

8. Add into the column 500  $\mu$ l of washing solution with ethanol. Centrifuge 10 000 - 13 000 rpm for 30 - 45 seconds. Discard flow-through. Repeat this step one more time.
9. After discarding the flow-through, place empty columns into same collection tubes and centrifuge 10 000 to 13 000 rpm for 30 - 45 seconds.
10. Place the columns in a new 1.5 ml tubes, add 50-100  $\mu$ l of the Elution buffer into the center of column filter (If 50  $\mu$ l of the elution buffer is used, the maximum concentration will be obtained; if 100  $\mu$ l of the Elution buffer is used, the maximum yield will be obtained), incubate 1 min at room temperature. Centrifuge 10 000 - 13 000 rpm for 1 minute.

**Also important!**When DNA concentration measuring on a cuvette, a DNA sample should be diluted 10-20 times with TE buffer pH 8.5 or with an elution buffer (do not use water), otherwise the measurements will be inaccurate.

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Issued by TÜV SÜD Management Service GmbH